




COMPANION ANIMAL DIAGNOSTIC REPORTS



GENERATED BY EVIDENCE-BASED ARTIFICIAL INTELLIGENCE






Feline Kidney Size


19:39:40

SampleOwner / SampleCat / 20-Mar-2025



Kidney size is assessed by forming the ratio of kidney length over the length of the L2 vertebra. This measurement of kidney size, when compared to published norms, is one way to help assess the state of the kidneys in the cat. Use this data in combination with all relevant clinical information.

Neutered



The right kidney is small with a kidney to L2 ratio of 1.90.


The left kidney is small with a kidney to L2 ratio of 1.81.

Small kidneys are consistent with chronic renal disease. Clinical significance should be correlated with bloodwork. A special diet could be considered:

Hill's Prescription Diet k/d

Reference: [Shiroma et al, 1999](#)

[Algorithm Notes](#)



Interpretation and Notes

If both kidneys are smaller than normal in your cat, it could suggest chronic kidney disease (CKD) which is a very common condition in older cats. However, radiographs only provide limited information, and further diagnostic steps are essential.



What These Results Might Indicate



Chronic Kidney Disease (CKD) is quite likely

- * Chronic degeneration and fibrosis cause the kidneys to shrink over time.
- * Kidneys may become: Smaller, Irregular in shape, Less distinct in margins
- * Common in older cats, especially over age 10.
- * Early CKD may still show normal kidney size, so small kidneys may imply advanced damage.



Other (less common) causes of bilaterally small kidneys:

- * Congenital hypoplasia (kidneys never developed fully; rare)
- * End-stage renal disease from previous toxin or infection
- * Post-inflammatory scarring (e.g. from feline infectious peritonitis)



Recommended Follow-Up



Bloodwork

- * BUN, Creatinine, and SDMA to assess global kidney function
- * Check electrolytes, phosphorus, potassium, etc.



Urinalysis

- * Look for concentrating ability (USG), protein, or signs of infection
- * UPC ratio if proteinuria is found



Blood Pressure Check

- * Hypertension is common in cats with kidney changes



Abdominal Ultrasound (Highly Recommended)

- * Gives a detailed look at internal kidney structure
- * May provide evidence of hydronephrosis, stones, tumors, or other issues



Prognosis / Action



Lab work rules out CKD → Explore less common causes, including congenital issues




Other tests confirm CKD → Chronic and irreversible, but manageable, especially in early stages. Special renal diet is recommended.



Feline Kidney Size


18:39:40

SampleOwner / Kitty / 26-Jan-2025



Kidney size is assessed by forming the ratio of kidney length over the length of the L2 vertebra. This measurement of kidney size, when compared to published norms, is one way to help assess the state of the kidneys in the cat. Use this data in combination with all relevant clinical information.

Unknown Neutered Status



The right kidney size is normal with a kidney to L2 ratio of 2.10.


The left kidney is small with a kidney to L2 ratio of 1.88.

Small kidneys are consistent with chronic renal disease. Clinical significance should be correlated with bloodwork. A special diet could be considered:

Hill's Prescription Diet k/d

Reference: [Shiroma et al, 1999](#)

[Algorithm Notes](#)







Interpretation and Notes

If your cat has one kidney that is smaller than normal and the other kidney is normal in size, it raises several possible interpretations. This asymmetry in kidney size is a clinically important finding and can point to a variety of underlying conditions.







What These Results Might Indicate

-  Chronic Kidney Disease (CKD) affecting one kidney. The smaller kidney may be scarred or atrophied due to:
 - * Past infection, inflammation, or trauma
 - * Obstruction (like a ureteral stone) that caused damage over time
 - * The normal-sized kidney may be compensating, doing more of the filtration work
 - * Common in middle-aged to older cats
-  Congenital Hypoplasia (underdeveloped kidney). The small kidney may have never developed fully.
 - * The other kidney typically enlarges (compensatory hypertrophy) to maintain function
 - * Often asymptomatic unless the functioning kidney becomes diseased later
-  Previous Kidney Injury. From ischemia (poor blood supply), toxins, or trauma.
 - * Causes permanent shrinkage and loss of function in one kidney
-  Chronic Obstructive Uropathy. Past or current ureteral obstruction (stone, stricture) can cause long-term damage
 - * A history of straining to urinate, vomiting, or pain may be present







Recommended Follow-Up

-  Bloodwork
 - * BUN, Creatinine, and SDMA to assess global kidney function
 - * Check electrolytes, phosphorus, potassium, etc.
-  Urinalysis
 - * Look for concentrating ability (USG), protein, or signs of infection
 - * UPC ratio if proteinuria is found
-  Blood Pressure Check
 - * Hypertension is common in cats with kidney changes
-  Abdominal Ultrasound (Highly Recommended)
 - * Gives a detailed look at internal kidney structure
 - * May provide evidence of hydronephrosis, stones, tumors, or other issues



Prognosis / Action



-  Small kidney + normal lab work  Likely past injury or congenital issue. Recheck periodically.
-  Small kidney + abnormal bloodwork  Likely unilateral CKD with possible early disease. Suggest a special renal diet.



Canine Abdomen

11:47:12

Smith / Yodo / 8-Jul-2022



Owner: PATTERSON, BONNIE
 Animal: MAGGIE
 Patient Age: 10
 Date: 18-Apr-2019 11:45:46
 R Abdomen Lateral
 Crater Road Veterinary

T13-L1	<input checked="" type="checkbox"/>	Pregnant	<input checked="" type="checkbox"/>	Urethral Calculi	<input checked="" type="checkbox"/>
L1-L2	<input checked="" type="checkbox"/>	Barium Swallow	<input checked="" type="checkbox"/>	Visible Prostate	<input checked="" type="checkbox"/>
L2-L3	<input checked="" type="checkbox"/>	Mild GD	<input checked="" type="checkbox"/>	Enlarged Colon	<input checked="" type="checkbox"/>
L3-L4	<input checked="" type="checkbox"/>	GD or GDV	<input checked="" type="checkbox"/>	Rotation	<input checked="" type="checkbox"/>
L4-L5	<input checked="" type="checkbox"/>	Enlarged Liver	<input checked="" type="checkbox"/>	(OK= <input checked="" type="checkbox"/> Found= <input checked="" type="checkbox"/> N/A= <input type="checkbox"/> Info= <input type="checkbox"/>)	
L5-L6	<input checked="" type="checkbox"/>	Micro Liver	<input checked="" type="checkbox"/>		
L6-L7	<input checked="" type="checkbox"/>	Renal Mineralization	<input checked="" type="checkbox"/>		
L7-S1	<input type="checkbox"/>	Bladder Calculi	<input checked="" type="checkbox"/>		

Algorithm Notes


References: [J. Choi et al., 2013](#)



Cardiac Thoracic Ratio


11:47:12

Smith / Yodo / 8-Jul-2022



The Cardiac Thoracic Ratio measurement is one way to help assess the size of the canine heart. It measures the size of the heart as a fraction of the width of the thorax at the 8-th rib. Use this result in combination with all relevant clinical information. Patient rotation less than 5% is best.

Cardiac Thoracic Ratio (CTR) for 14,231 Dogs




The CTR measurement for this dog is 0.501. This value is indicated on the histogram to the left by the vertical black line. Large values might indicate cardiomegaly, small values might indicate microcardia.


The CTR is within the central 70% of dogs.

Reference: [Azevedo et al, 2016](#)

Algorithm Notes



Elbow Health Score: 74%




References for understanding these measures

HRI: [A. Mostafa et al., 2019](#)

STS: [T.J. Smith et al., 2009](#)


Several measurements have been made on this image. You can learn more about each of these measures by clicking on the blue links in this report. In order to give an overall numeric assessment, we have computed the 'Elbow Health Score' based on what is visible in this radiograph. This measure, serves as one data point for the veterinarian, to help monitor the pet's well-being.




The 'Elbow Health Score' for this elbow is 74%. This value is indicated on the histogram to the left by the vertical black line. The further to the right on this diagram, the better.

This elbow is in the 3rd (below ave) quartile.

[Algorithm Notes](#)



The Norberg Hip Angle measurements are one way to help assess the morphology of canine hip joints. Each hip joint is assigned a score - the larger the value, the better. The measurement is used as an indicator in combination with other relevant clinical information.



The Norberg measurements for this pet are 109.2 for the right hip and 110.9 for the left hip. These values are indicated on the histogram to the left by the vertical black lines. The further to the right on this diagram, the better.

The right hip is in the 1st (best) quartile.

The left hip is in the 1st (best) quartile.

Reference: [M. Flückiger, 2007](#)

[Algorithm Notes](#)


Distribution based on 29,000 canine hips




Canine Shoulder

10:13:18

Smith / Yodo / 8-Jul-2022



Original Image



Zoom in on Shoulder

The area around the shoulder joint has been checked for various anomalies. If an anomaly was found, it is circled in red on the large image above, and is highlighted in red on the indicator below. Green indicators below signify that the named anomaly was not found. Yellow indicators signify that the determination could not be made for this image.


Crater Lesion	■	Osteochondroma	■
Caudal Capsule Osteophyte	■	Humeral Lysis	■
Caudal Sac Osteophyte	■	Frac Supraglenoid Tubercle	■
Caudal Capsule Ossicle	■	Humeral Sclerosis	■
Caudal Sac Ossicle	■	(OK= ● Found= ● N/A= ● Info= ●)	
Bicipital Tendon Anomaly	■		



Canine TPLO Report


11:47:12

Smith / Yodo / 8-Jul-2022



The Tibial Plane Leveling Osteotomy (TPLO) is a procedure to lower the effective value of the Tibial Plane Angle (TPA). High values of TPA are weakly correlated with risk of CCL damage. Other stifle procedures might also be considered.

TPA



The TPA measurement for this stifle is 11.8 degrees and is indicated on the diagram to the left by the blue arrow.

The Tibial Plane Angle (TPA) of the stifle joint is normal or low.

For best results, position the leg with the stifle flexed so that 'Leg Angle' is 90 degrees.

To calibrate for saw sizing, please use our ['AutoScaler'](#)

Ref: [Todorovic et al, 2022](#)


[Algorithm Notes](#)




Canine TTA Report

11:47:12

Smith / Yodo / 8-Jul-2022



The Tibia Tuberosity Advancement (TTA) is one measurement made to help assess whether this stifle might benefit from the TTA procedure. If the measured advancement is large, the TTA procedure might be of benefit. Other stifle procedures might also be considered.



The TTA measurement for this stifle is 3.9 degrees and is indicated on the histogram to the left by the vertical black line. The green-colored region of the graph represents the 'normal' range.

The relative positioning of the bones of the stifle joint appear normal. This stifle probably would not benefit from the TTA procedure.

Ref: [Background Info on TTA](#)

Algorithm Notes


Distribution for dogs under veterinary care




Canine VHS/VLAS Report

11:47:12

Smith / Yodo / 8-Jul-2022



The Vertebral Heart Size (VHS) and Vertebral Left Atrial Size (VLAS) are two objective measurements of heart size, which in combination with other clinical data, can help stage acquired canine heart disease including myxomatous mitral valve disease (MMVD) and dilated cardiomyopathy (DCM).



In the absence of an echocardiogram, older small breed dogs with moderately loud murmurs consistent with mitral regurgitation can be considered Stage B2 if the VLAS is ≥ 3.0 and or the VHS is ≥ 11.5 . In large breed dogs suspected to have heart disease, VHS ≥ 11.5 is strongly supportive of clinically relevant cardiomegaly and warrants additional evaluation.

References for understanding these measures

Ref: [The Cardiac Education Group](#)

Ref: [The ACVIM MMVD guidelines](#)

Ref: [Large Breed Dogs](#)

Breed-Specific for norwich terrier

VHS Normal Range: 9.4 - 11.8 (Mean +/- 2 SD)

Ref: [Colin Taylor et al \(2020\)](#)

Analysis powered by www.MetronMind.com

[Algorithm Notes](#)

These results should be considered in light of all relevant clinical information.