


# Excellence in Digital Radiography. EMPOWER1717

The exclusive features and functions of the Metron DVM software set the Empower 1717 apart from other systems in the industry. Amazing resolution and dynamic range, paired with the largest digital radiography panel offered in Veterinary medicine today, insures a complete image on even the largest patient. Includes full 3-year warranty, 3 copies of the Metron DVM software and one year of Internet-based remote storage from \*Asteris.biz.

- Instantly capture images in less than 4 seconds
- Create DICOM Compliant CDs and DVDs
- Store up to 55,000 images
- Easy to use database for retrieval and organization
- 24 hour Help Desk



Reclaim Valuable Technician's Time



Evaluate X-Rays in Seconds



Detailed Images



Panel Mounts into X-Ray Table

## Metron DVM Software

- ▶ Designed exclusively for veterinary use
- ▶ Intuitive, easy-to-learn features for doctors and staff
- ▶ PACS feature is standard, not extra
- ▶ Generates reports and DICOM-compliant files
- ▶ Special features include heart score, TPLO, guided mark up, dental templates, relative bone density measurement and magnification corrected anatomic measurement tools

## Technical Specifications

Panel Information	17" x 17" Empower panel
Scintillator	CSI
Contrast Resolution	14 bit (16,384 shades of gray)
Spatial Resolution	143 micron - 3.5 lp/mm
Image Server	1 TB RAID 1 - 4 GB RAM
Image Storage	Up to 55,000 images
Software License	3 copies of Metron DVM Software
Training	1 day on-site training included
Warranty Information	3 year incl. software updates & upgrades

Experience the ease and performance of the Empower 1717 by scheduling a hands-on demonstration in your office. Our experts will even run Return on Investment figures to show you how soon you'll be saving money and time. Call today!



Best Quality. Best Value. Best of Breed.

**800.458.8890**

[www.veteldiagnostics.com](http://www.veteldiagnostics.com)