

# Excellence in Digital Radiography.

# EMPOWERWIRELESS

The New WIRELESS panel on this companion animal system offers the highest quality radiographs available on the market today, while allowing unimagined flexibility in capturing images. Gain the ability to acquire radiographic images anywhere in your hospital, saving time and patient movement. This system's portability, user-friendliness and quality are designed with maximum efficiency and value. It includes a full 3-year warranty, 1 copy of the Vetel InSight™ software, and one full year of internet-based image storage from \*Asteris.biz. Empower your practice with confidence and excellence.

**NEW!**



Wireless Panel with Charging Base

- **WIRELESS PANEL!**  
14" X 17" imaging area
- 64,000 shades of gray
- Instantly capture images in less than 9 seconds
- Create DICOM compliant CDs and DVDs
- Store up to 55,000 images
- Easy to use database for retrieval and organization
- 24 hour Help Desk



\*Cart Optional



Reclaim Valuable Technician's Time



Evaluate X-Rays in Seconds



Detailed Images

## Vetel InSight™ Software

- ▶ Multi-planar bone density analysis – analyze up to four areas at one time
- ▶ Image processing algorithms that maximize spatial resolution and dynamic range regardless of x-ray generator milliamperage
- ▶ Exclusive monitor calibration capabilities
- ▶ Voice prompting for proper sequencing of shots
- ▶ Email images directly from InSight's address book
- ▶ Ability to generate reports and DICOM-compliant files
- ▶ Capacity to customize records and attach images to important client and patient information

## Technical Specifications

Panel Information	14" x 17" Empower wireless panel
Scintillator	CSi
Contrast Resolution	16 bit (64,000 shades of gray)
Spatial Resolution	144 micron - 3.5 lp/mm
Image Server	1 TB RAID 1 - 4 GB RAM
Image Storage	Up to 55,000 images
Software License	1 copy of Vetel InSight™ Software
Training	1 day on-site training included
Warranty Information	3 year incl. software updates & upgrades

Experience the ease and performance of the Empower Wireless by scheduling a hands-on demonstration in your office. Our experts will even run Return on Investment figures to show you how soon you'll be saving money and time. Call today!



Best Quality. Best Value. Best of Breed.

**800.458.8890**  
www.veteldiagnostics.com