



Excellence in Digital Radiography.

EXPLORERCA

The most popular digital radiography panel is available in an easy to use companion animal system. The affordable ExplorerCA combines high performance with sturdy reliability. Powerful Metron DVM software provides exclusive modules and diagnostic features to complement Vetel Diagnostic's leading DR panel manufactured by Canon®.

- Instantly capture images in less than 4 seconds
- Store up to 23,000 images
- Easy to use database for retrieval and organization
- 24 hour Help Desk



Reclaim Valuable Technician's Time



Evaluate X-Rays in Seconds



Detailed Images



Easy to Use Portable Panel

Metron DVM Software

- ▶ Designed exclusively for veterinary use
- ▶ Intuitive, easy-to-learn features for doctors and staff
- ▶ PACS feature is standard, not extra
- ▶ Generates reports and DICOM-compliant files
- ▶ Special features include heart score, TPLO, guided mark up, dental templates, relative bone density measurement and magnification corrected anatomic measurement tools

Technical Specifications

Panel Information	14" X 17" Canon® 50G panel
Scintillator	GadOX
Contrast Resolution	12 bits (4,096 shades of gray)
Spatial Resolution	160 micron - 3.1 lp/mm
Image Server	320 GB RAID 1 - 2GB RAM
Image Storage	Up to 23,000 images
Software License	1 copy of Metron DVM Software
Training	1 day on-site training included
Warranty Information	1 year incl. software updates & upgrades

Experience the ease and performance of the ExplorerCA by scheduling a hands-on demonstration in your office. Our experts will even run Return on Investment figures to show you how soon you'll be saving money and time. Call today!



Best Quality. Best Value. Best of Breed.

800.458.8890
www.veteldiagnostics.com