



Excellence in Digital Radiography. EMPOWEREQ

A demanding on-the-go practice needs a quality, compact solution. With the New WIRELESS panel—offering the highest quality images available on the market today—you will experience unimagined flexibility and freedom from troublesome cables. This user-friendly system includes all essential equine accessories, 3-year warranty, 1 copy of the Vetel InSight™ software and 1 full year of internet-based image storage from *Asteris.biz. The EmpowerEQ combines unsurpassed image quality, portability and ease of use. Hit the road with confidence and excellence.

- **WIRELESS 14" x17" Panel!**
- Instantly capture images in less than 9 seconds
- 64,000 shades of gray
- Store up to 50,000 images
- Easy to use database for retrieval and organization
- 24 hour Help Desk



Rugged and Portable



Easy to Use



Wireless Docking Station Available



Wireless Panel Fits Securely Inside

Vetel InSight™ Software

- ▶ Multi-planar bone density analysis – analyze up to four areas at one time
- ▶ Image processing algorithms that maximize spatial resolution and dynamic range regardless of x-ray generator milliamperage
- ▶ Exclusive monitor calibration capabilities
- ▶ Voice prompting for proper sequencing of shots
- ▶ Email images directly from InSight's address book
- ▶ Ability to generate reports and DICOM-compliant files
- ▶ Capacity to customize records and attach images to important client and patient information

Technical Specifications

System Weight	35 lbs
System Size	11" x 22" x 19"
Hard Drive	500 GB RAID 1
RAM	4 GB
Image Storage	Up to 50,000 images
Software Licenses	1 copy of Vetel InSight™ Software
Training	1 day on-site training
Warranty Information	3 year including software upgrades
Panel	16 bit - 14" x 17" imaging area

Experience the ease and performance of the EmpowerEQ by scheduling a hands-on demonstration in your office. Our experts will even run Return on Investment figures to show you how soon you'll be saving money and time. Call today!



Best Quality. Best Value. Best of Breed.
800.458.8890
www.veteldiagnostics.com