



# Excellence in Digital Radiography.


## ENVISION1210

This powerful equine package has everything you need for maximum value and efficiency. The full-featured Envision1210 digital radiography system includes a collection of indispensable accessories, 3-year warranty, 3 copies of Metron DVM software and one full year of internet-based image storage. High performance system with 17" display, RAID-1 hard drive and 4GB of RAM. The Envision1210 is the out-of-the-box solution offering fantastic value and features for any busy equine practice.

- Samsung® CSI Panel
- Instantly capture images in less than 4 seconds
- 16,384 shades of gray
- Store up to 50,000 images
- 10" x 12" imaging area
- 24 hour Help Desk



Rugged and Portable



Easy to Use



Durable Accessories



Panel Fits Securely Inside Case

### Metron DVM Software

- ▶ Designed exclusively for veterinary use
- ▶ Intuitive, easy-to-learn features for doctors and staff
- ▶ Image processing algorithms maximize dynamic range of images from x-ray generators, ultrasound and endoscopy
- ▶ PACS feature is standard, not extra
- ▶ Generates reports and DICOM-compliant files
- ▶ Special features include application specific modules, relative bone density measurement, 3-D rendering, and time-lapse analysis morphing

### Technical Specifications

System Weight	21 lbs
System Size	21" x 14" x 7"
Hard Drive	500 GB RAID 1
RAM	4 GB
Image Storage	Up to 50,000 images
Software Licenses	3 copies of Metron DVM Software
Training	1 day on-site training
Warranty Information	1 year including software upgrades
Panel	14 bit - 10" x 12" imaging area

Experience the ease and performance of the Envision1210 by scheduling a hands-on demonstration in your office. Our experts will even run Return on Investment figures to show you how soon you'll be saving money and time. Call today!



Best Quality. Best Value. Best of Breed.

**800.458.8890**

[www.veteldiagnostics.com](http://www.veteldiagnostics.com)