

Excellence in On-Site Thermal Imaging.

THERMACAMB360


Exceptional resolution and FLIR's exclusive Advanced Signal Processing yield razor-sharp image quality. The B360 packs all your "wish list" items into one powerful unit with features and functions to make identification and interpretation of infrared images even easier. Exceptional optics, large 3.5" display, long battery life and advanced analysis tools are just a few of the B360's distinctive elements.




- 3.5" Color Touch-screen LCD
- Thermal Fusion Functionality
- Removable SD/Memory Card, USB & Video Out
- Independently Rotatable Optics Head & Display Screen
- Built-in, patented Laser LocatIR™
- Hot/Cold Spot & Visual Alarms
- JPEG Format, Email-friendly




Intuitive Interface



Easy to Use



Hoof Imbalance



Shoulder Inflammation

FLIR Infrared Cameras

- ▶ FLIR infrared cameras incorporate 30+ years of experience in thermal imaging.
- ▶ You get superior ergonomics, advanced signal processing, ambient temperature control and a host of infrared camera and thermal imaging exclusives.
- ▶ FLIR offers the widest selection of IR cameras and thermal imagers, the best post-sale technical support, the world's largest and most experienced thermal imaging training organization—the Infrared Training Center—and award-winning customer service.

Technical Specifications

Camera Weight	1.94 lbs (0.88kg)
Camera Size	4.2" x 7.9" x 4.9"
Field of View / Min. Focus	25° x 19° / 0.4m (1.31ft)
Thermal Sensitivity	0.06°C @ + 30° (+86°F) / 60mK
Spectral Range	7.5 to 13µm
Image Storage	1000+ JPEG images
Training	Online Support Available
Warranty Information	Various Plans Available

Join us for our semi-annual, 4 day, CE approved training seminars. Visit our web site for dates and details: www.veteldiagnostics.com

Experience the ease and performance of the ThermoCam B360. Our experts will even run Return on Investment figures to show you how soon you'll be saving money and time. Call today!



Recommended by a veterinarian, for veterinarians.

800.458.8890
www.veteldiagnostics.com