



2024 EQUIPMENT CATALOG

PREMIER DIAGNOSTIC SOLUTIONS FOR VETERINARIANS

HIGH DEFINITION DIGITAL RADIOGRAPHY: PG. 4

ULTRASOUND SOLUTIONS FOR YOUR CLINIC: PG.8

TABLE OF CONTENTS

PRACTICE WISH LIST	3
DIGITAL RADIOGRAPHY	4
DIGITAL DENTAL RADIOGRAPHY	6
ULTRASOUND - DRAMINSKI	8
PALM-SIZED ULTRASOUND - ECHO IQ	10
ULTRASOUND - CHISON	11
THERMOGRAPHY	12
X-RAY GENERATORS AND TABLES	14
NEEDLE ARTHROSCOPY	15
ENDOSCOPY	16
AI ASSISTANT AND REPORT GENERATOR	17
ONLINE EDUCATION	18
WHY VETERINARIANS LOVE WORKING WITH VETEL	19
FINANCING	20

SERVING THE VETERINARY COMMUNITY FOR OVER 25 YEARS



WWW.VETELDIAGNOSTICS.COM

800-458-8890 EXT. 3

INFO@VETELDIAGNOSTICS.COM

VETEL DIAGNOSTICS IS A TRUSTED PARTNER FOR VETERINARY SUCCESS AND PATIENT WELLNESS ACROSS NORTH AMERICA. VETEL WORKS WITH VETERINARY EXPERTS TO CUSTOMIZE IMAGING EQUIPMENT FOR ALL SPECIALTIES. IMPROVE YOUR PATIENT CARE WITH STATE-OF-THE-ART DIAGNOSTIC MODALITIES, ENABLING EFFICIENT AND ACCURATE RESULTS EVERY TIME. DESIGNED FOR VETERINARIANS BY VETERINARIANS, VETEL EQUIPMENT STANDS AT THE TOP OF ITS CLASS SIMPLY BY MEETING YOU WHERE YOU NEED US.

PRACTICE WISH LIST

DIGITAL RADIOGRAPHY

ULTRASOUND

THERMAL CAMERAS

DENTAL DIGITAL RADIOGRAPHY


ENDOSCOPY

X-RAY GENERATORS AND TABLES

ARTHROSCOPY

NEED SOMETHING ELSE? WE'LL HELP FIND YOUR SOLUTION.

EMPOWER HDR™



COMPANION ANIMAL DIGITAL RADIOGRAPHY EXCELLENCE FOR EVERY PATIENT

- CLEAR HIGH DYNAMIC RANGE IMAGING FOR PRECISE AND ACCURATE DIAGNOSIS
- POWERFUL DELL ALL-IN-ONE PC
- GENERATE AUTOMATED MEASUREMENTS AND REPORTS
- EXTRA PROTECTION WITH MISHAP IMMUNITY
- COMPLETE ACQUISITION AND PACS SOFTWARE SUITE
- LUMINA IQ™ AI-POWERED INTELLIGENT IMAGE PROCESSING
- 14" X 17" AND 17" X 17" PANEL SIZES COMPATIBLE WITH STANDARD BUCKY TRAYS



FIVE YEAR WARRANTY

SOFTWARE AND HARDWARE PROTECTION PLAN

[Click here to view this product online](#)

EMPOWER HD™



PREMIUM COMPANION ANIMAL DIGITAL RADIOGRAPHY ADVANCED IMAGING TO SUPPORT YOUR DIAGNOSIS

- HIGH DEFINITION TECHNOLOGY FOR A CRYSTAL CLEAR IMAGE
- FULL-FIELD AED PANELS
- MANAGE PATIENT STUDIES SEAMLESSLY WITH A POWERFUL DELL ALL-IN-ONE PC
- FAST 1.5 SECOND IMAGE CAPTURE AND PROCESSING
- COMPLETE ACQUISITION AND PACS SOFTWARE SUITE
- LUMINA IQ™ AI-POWERED INTELLIGENT IMAGE PROCESSING
- 14" X 17" AND 17" X 17" PANEL SIZES COMPATIBLE WITH STANDARD BUCKY TRAYS



FIVE YEAR WARRANTY

SOFTWARE AND HARDWARE PROTECTION PLAN

ENVISION DENTAL™

DIGITAL DENTAL RADIOGRAPHY



CUSTOMIZABLE DENTAL SUITES FOR CLINICS AND MOBILE PRACTICES

- WALL MOUNTED X-RAY GENERATORS
- HANDHELD X-RAY GENERATORS
- CHOICE OF DENTAL SENSOR SIZE #1, #2, OR #3
- ADAPTABLE DELL 2-IN-1 PC
- COMPREHENSIVE DENTAL ACQUISITION SOFTWARE INCLUDING 3D IMAGE AND POSITIONING GUIDE



DENTAL IMAGING SIMPLIFIED

[Click here to view this product online](#)

ENVISION EXL[®] DENTAL CR


MIXED PRACTICE DIGITAL DENTAL IMAGING




MINIMUM SIZE, MAXIMUM PERFORMANCE.

- DENTAL CR SYSTEM WITH ALL THE FLEXIBILITY AND COMFORT OF FILM WITH THE IMMEDIACY OF DIGITAL RADIOGRAPHY
- SMART LOAD-TOP DESIGN HARNESSES GRAVITY FOR SIMPLE AND FUSS-FREE USE
- TINY ENOUGH TO FIT INTO CROWDED DENTAL OPERATORIES
- METRON ACQUISITION SOFTWARE INCLUDING 3D IMAGE AND POSITIONING GUIDE
- EXCLUSIVE CLEANING MODE FOR EASY MAINTENANCE

PATIENT-FIRST DENTAL IMAGING.



BLUE DRAMIŃSKI ULTRASOUND



A PREMIER PORTABLE ULTRASOUND FOR COMPREHENSIVE EXAMS

- HIGHLY RESPONSIVE INTUITIVE INTERFACE
- VETERINARIAN APPROVED PRESETS
- IDEAL FOR 24/7 HOSPITALS AND PRACTICES
- SUPPORTS ANY VETERINARY SPECIALTY
- COLOR, POWER, AND PULSE WAVE DOPPLER



CONTRAST AND
SHARPNESS
ENHANCEMENT



GREYSCALE CURVE
MATCHED TO THE
HUMAN EYE



TWO YEAR WARRANTY

DRAMINSKI SYSTEMS. ONE YEAR ON PROBES.

[Click here to view this product online](#)



WIDER FIELD
OF VIEW

iSCAN 2 MULTI

DRAMIŃSKI ULTRASOUND



- MULTI-TASK IN CLINIC OR MOBILE SETTINGS WITH A WIDE SELECTION OF PROBES
- UNDER 6 LBS WITH BATTERY
- WORK CONTINUOUSLY FOR UP TO 7 HOURS ON A FULL CHARGE
- CLEAR 7" DISPLAY
- PROBES INTERCHANGEABLE BETWEEN iSCAN 2 MULTI AND BLUE

DRAMIŃSKI PROBE SELECTION



CONVEX ABDOMINAL



LINEAR ABD / MSK



MICROCONVEX ABDOMINAL



CONVEX RECTAL



LINEAR RECTAL



MICROCONVEX RECTAL




TWO YEAR WARRANTY

DRAMINSKI SYSTEMS. ONE YEAR ON PROBES

[Click here to view this product online](#)

ECHO IQ WIRELESS ULTRASOUND

[Click here to view this product online](#)



A POWERFUL PALM-SIZED ULTRASOUND ALWAYS WITHIN YOUR REACH

POWER IN SMALL PACKAGES

- WIRELESS 5GHz CONNECTION
- RELIABLE HIGH-DEFINITION IMAGING
- COMFORTABLY FITS IN YOUR POCKET

ONE DEVICE, ENDLESS APPLICATIONS

- HOT-SWAP PROBES MID-EXAM WITH 1-CLICK
- CHOICE OF PROBES AVAILABLE FOR ANY SPECIALTY



*CONVEX, MICROCONVEX,
LINEAR AND RECTAL PROBE
TYPES AVAILABLE*

SONOBOOK[®] ULTRASOUND

[Click here to view this product online](#)



SONOBOOK 6

- COMPACT LAPTOP DESIGN FOR EASY MOBILITY
- LESS THAN 12 LBS
- CRISP IMAGING WITH A 15" HIGH-RES LED DISPLAY
- SPECIES-SPECIFIC PRESETS
- USER FRIENDLY INTERFACE WITH SHORTCUT KEY FOR QUICK ACCESS

SONOBOOK 8 VET

- PREMIUM CLASS ULTRASOUND WITH A 15" LED HIGH-RES DISPLAY
- COMPREHENSIVE CARDIAC MEASUREMENT PACKAGES
- FHI HARMONIC IMAGING
- AUTOMATIC MEASUREMENT AND IMAGE OPTIMIZATION
- SUPERNEEDLE ENHANCEMENT
- 12 PROBE OPTIONS
- ULTRA-REMOTE LIVE WEB SHARING



THREE YEAR WARRANTY
CHISON ULTRASOUND SYSTEMS



ULTRA-REMOTE LIVE WEB SHARING


REMOTE INSTRUCTION AND CONSULTING MADE EASY.
COLLABORATE WITH YOUR COLLEAGUES IN REAL-TIME.

[Click here to view this product online](#)



THERMOGRAPHY

THE VETERINARIANS SIXTH SENSE



OFFER THERMAL SERVICES & GROW YOUR PRACTICE

- INCREASE CLIENT CONFIDENCE BY SHOWING REAL-TIME THERMAL ARTIFACTS AND IMPROVEMENTS IN MUSCULOSKELETAL HEALTH OVER TIME
- QUICKLY LOCATE INFLAMMATION, PAINFUL AREAS, OR NERVE ROOT ARTIFACTS
- PRACTICE PREVENTATIVE CARE BY IDENTIFYING NON-SYMPTOMATIC INDICATORS OF LAMENESS EARLY ON

REFRESH YOUR THERMOGRAPHY KNOWLEDGE

STUDY VETERINARY THERMOGRAPHY ONLINE AT YOUR OWN PACE WITH VETEL'S COMPREHENSIVE EDUCATION COURSE TAUGHT BY BOARD-CERTIFIED EXPERTS IN THERMOGRAPHY

THERMAL CAMERAS

PREMIUM CLASS THERMAL IMAGING ONLY FROM



- T-530**
- 320 X 240 IR RESOLUTION
 - THERMAL SENSITIVITY: -20°C - 650°C
 - 4" LCD TOUCH DISPLAY WITH 180° ROTATING LENS
 - VIDEO RECORDING AND STREAMING

- T-560**
- 640 X 480 IR RESOLUTION
 - THERMAL SENSITIVITY: -20°C - 1500°C
 - 4" LCD TOUCH DISPLAY WITH 180° ROTATING LENS
 - PICTURE-IN-PICTURE CAPABILITIES



ALL MODELS WIFI-CAPABLE
AND PLUG INTO ANY PC FOR
EASY DATA MANAGEMENT

[Click here to view these products online](#)

E76

- 320 X 240 IR RESOLUTION
- THERMAL SENSITIVITY: <40 MK AT 30°C
- 4" LCD TOUCH DISPLAY WITH AUTO-ROTATION
- LASER-ASSISTED GUIDANCE AND AUTO-FOCUS
- VOICE ANNOTATION



FLIR TOOLS WITH ULTRAMAX®

FLIR TOOLS IS A FREE SOFTWARE THAT ALLOWS THERMAL CAMERA USERS TO QUICKLY IMPORT, EDIT, AND ANALYZE IMAGES, TURNING THEM INTO PROFESSIONAL REPORTS. INCLUDED IS ULTRAMAX IMAGE ENHANCEMENT, WHICH MULTIPLIES PIXELS IN A TARGET AREA FOR SUPERIOR IMAGE DETAIL AND DEFINITION.

X-RAY GENERATORS

AV CHOICE HIGH FREQUENCY GENERATOR

SINGLE OR THREE - PHASE



- 30KW, HIGH FREQUENCY, SINGLE PHASE, 400MA, 125KVP, WITH ANATOMICAL PROGRAMMING CONSOLE
- OPTIONAL 32KW, HIGH FREQUENCY, THREE PHASE, 600MA, 125 KVP, WITH ANATOMICAL PROGRAMMING CONSOLE AVAILABLE
- AVAILABLE WITH STATIONARY, 2-WAY FLOAT, OR 4-WAY FLOAT TOP TABLE
- URINE TRAPS FOR EASY CLEAN UP

EV-800 ELEVATED FLOAT TOP TABLE

- SUPPORTS UP TO 800 POUNDS
- 36" WIDE COMPLETELY FLAT TABLE TOP
- 86" LENGTH
- INTEGRAL COLLISION SENSORS
- ELEVATING RANGE 22.25" TO 34.25"
- STANDARD TABLE-SIDE CONTROL
- OPTIONAL CARBON/FIBER TOP



UMG/DEL MEDICAL

[Click here to view these products online](#)

ARTHROSCOPY

NEEDLEVIEW SYSTEM

- ARTHROSCOPE SUITE WITH A VISUALIZATION TOOL THE SIZE OF AN 18 GAUGE NEEDLE
- PERFORM DIAGNOSTIC AND THERAPEUTIC PROCEDURES
- COMPACT, PORTABLE, AND FULLY DIGITAL

SMALL ANIMAL APPLICATIONS

ELBOW, HIP, STIFLE, SHOULDER



CANINE JOINT INJECTION

ASK ABOUT TRAINING



TELEVIEW USB ENDOSCOPES



*COMPUTER AND TRIPOD NOT INCLUDED

- **HIGH DEFINITION CAMERA** WITH SUPER BRIGHT LED ILLUMINATION TO GUIDE YOUR WAY
- IRRIGATION CHANNEL KEEPS VIEW UNOBSTRUCTED
- WORKING CHANNEL FOR BIOPSIES AND FOREIGN BODY REMOVAL
- AIR-WATER CHANNEL FOR LAVAGE OR INSUFFLATION



LENGTHS: 1M, 1.5M, 3M, 3.5M & CUSTOM

DIAMETERS: 8.7MM, 9.8MM, & 12.6MM

[Click here to view this products online](#)



H1 MODEL



- THREE WAVELENGTHS - 650NM, 810NM, 915NM
- 15 WATTS - HIGH POWER ALLOWS DEEPER PENETRATION OF TISSUE IN LESS TIME
- SOLID STATE LASER DIODE MODULE
- DEVICE WEIGHT - 4.4 LBS
- 7 INCH 1920X1200 HD TOUCH SCREEN

LIFETIME DIODE WARRANTY - 3 YEAR MANUFACTURERS WARRANTY




PLATINUM P4

- FOUR WAVELENGTHS - 650NM, 810NM, 980NM
- 24 WATTS - MORE POWERFUL LASER FOR SHORTER TREATMENTS AND FASTER RESULTS
- WIFI PERFORMANCE TRACKING APP
- OVER 1300 EXPERT TREATMENT PROTOCOLS
- 7 INCH 720P HD TOUCHSCREEN
- SOLID STATE LASER DIODE MODULE

NEED AN ASSISTANT?

EVIDENCE-BASED AI TECHNOLOGY STUDIES YOUR DIGITAL RADIOGRAPH AND
GENERATES COMPREHENSIVE DIAGNOSTIC REPORTS BASED ON NON-BIASED
ANALYTICS WITH JUST ONE CLICK



LET AI DO THE WORK FOR YOU

REVOLUTIONIZE YOUR WORKFLOW

RECLAIM LOST TIME

ASK US ABOUT METRON SOFTWARE

EDUCATION


VETEL HAS PUT TOGETHER ONLINE EDUCATION PROGRAMS FOR THERMOGRAPHY AND ULTRASONOGRAPHY WITH THE TOP EXPERTS IN EACH FIELD. WE PARTNER WITH MANUFACTURERS TO PROVIDE TRAINING FOR ALL THE TOOLS IN OUR DIAGNOSTIC TOOLBOX.

ONLINE COURSES INCLUDED FREE IN ALL ULTRASOUND OR THERMAL SYSTEMS

THERMOGRAPHY

- TH-01 - INTRODUCTION AND SAFETY AROUND ANIMALS
- TH-02 - ARTIFACTS, ENVIRONMENTAL CONCERNS & PHYSIOLOGY
- TH-03 - PERFORMING THE NORMAL EQUINE THERMOGRAPHIC EXAMINATION
- TH-04 - EQUINE LIMBS: THE FOOT, JOINTS, TENDONS AND LIGAMENTS
- TH-05 - BACK PAIN IN HORSES


+ 9 MORE MODULES



ULTRASOUND

- 01 - VETERINARY ULTRASOUND – INTRODUCTION AND PHYSICS
- 02 - VETERINARY ULTRASOUND – ARTIFACTS
- 03 - VETERINARY ULTRASOUND – HEPATOBILIARY SYSTEM
- 04 - VETERINARY ULTRASOUND – GALLBLADDER
- 05 - VETERINARY ULTRASOUND – SPLEEN

+ 6 MORE MODULES



[Click here to view EDU courses](#)

WHY VETERINARIANS LOVE WORKING WITH



"EXCELLENT PRODUCTS AND CUSTOMER SERVICE...THE SERVICE TEAM AMAZINGLY SUPPORTIVE AND KNOWLEDGEABLE!"

- JAN LUQUIRE, DVM

"VETEL UNDERSTOOD THE FRUSTRATIONS I HAD WITH MY OLD PORTABLE X-RAY UNIT. THEY WERE ABLE TO TAKE THAT UNIT AND SOFTWARE, FORMAT IT, AND SELL IT ON CONSIGNMENT. I AM HAPPY WITH MY TRANSITION TO VETEL'S RUGGED PORTABLE UNIT AND SOFTWARE."

- BUR CHURCH, DVM

"I HAVE BEEN A VETEL CUSTOMER FOR OVER 15 YEARS. QUALITY IMAGING EQUIPMENT WITH EXCELLENT 24/7 CUSTOMER SERVICE. THEY ABSOLUTELY UNDERSTAND THE NEEDS OF TODAY'S EQUINE PRACTITIONER."

- MARK AKIN, DVM

"VETEL WAS VERY HELPFUL HELPING ME TRANSITION FROM AN OLD INDIRECT SYSTEM TO THEIR SYSTEM EVEN THOUGH I AM ACROSS THE COUNTRY. IT WAS SEAMLESS AND THEY KEPT UP WITH ALL MY QUESTIONS! THEY ARE GREAT FOR SUPPORT AS WELL EVEN DOWN THE ROAD. . VERY HAPPY WITH THE SYSTEM. I HAVE A MOBILE PRACTICE FOR SMALL ANIMALS AND IT WORKS GREAT FOR OUR NEEDS."

- SARAH SKINNER, DVM

"CLARITY AND QUALITY OF IMAGE WAS BEST-IN-CLASS WITH THE RUGGED ENDURO UNIT, AND METRON SOFTWARE IS INTUITIVE, AND CUSTOMER SUPPORT... SECOND TO NONE! DON'T BE FOOLED BY THESE BIGGER COMPANIES. VETEL IS A GEM THAT MORE PEOPLE NEED TO DISCOVER. THANK YOU!"

- DRAFT COLT, DVM

"HAVING HAD OTHER COMPANIES XRAY SYSTEMS BEFORE I PURCHASED ONE FROM VETEL, I HAVE TO SAY VETEL SUPPORT IS FAR AND AWAY THE MOST RESPONSIVE AND EASIEST TO CONNECT WITH."

- HELEN POLAND, DVM

"WE HAVE LOVED WORKING WITH VETEL AND LOOK FORWARD TO MANY YEARS TO COME. COULD NOT SAY ENOUGH GOOD THINGS ABOUT THE STAFF, I TRULY FELT THAT THEY CARED AND WERE ON MY TEAM AS THEY WALKED THROUGH ALL THE PROBLEM SOLVING WITH ME. THANK YOU FOR HANGING IN THERE, AND GOING THE EXTRA MILE. WE ARE SO HAPPY WITH THE NEW MACHINE, AND SO GRATEFUL FOR ALL THE SUPPORT THAT CAME WITH IT ALONG THE WAY."

- BLACK DIAMOND EQUINE

FIRST YEAR FREE FINANCING

EXCLUSIVELY FROM



[Financing FAQ](#)

EXPERIENCE FREEDOM FROM FINANCING PAYMENTS
FOR THE FIRST 12 MONTHS OF YOUR TERM.*

HOW MUCH COULD YOU SAVE WITH
1 YEAR OF NO PAYMENTS?

800-458-8890

INFO@VETELDIAGNOSTICS.COM

*ON QUALIFYING EQUIPMENT PURCHASES OF \$10,000 OR MORE